

Transmitters

ZS-910/920/930/940PG


ZM-920/930/940PG

ZS-530PG


**Better care
by wireless**




Emergency



Nurse Station



Transport



Ward

Fighting Disease with Electronics

NIHON KOHDEN

High quality patient care by wireless

— patient is free from the bed —

❖ Wireless monitoring system

You can create a wireless monitoring system with transmitters and a central monitor with multiple patient receiver or a NIHON KOHDEN telemetry system.

❖ Increase patient safety

By continuous measurement of important vital parameters, patient condition is always monitored, quick reactions supported.

❖ Increased patient comfort

No cable between the patient and monitor, so the patient is free from the bed. The transmitter is compact and light-weight so it is easy to wear.

❖ More efficient patient care

Instead of visit and manual measurement, automatic measurement saves the care staff's time.

❖ Easy to use

No complicated preparation to start measurement. The CALL button activates the recorder and nurse call at the central monitor (ZS-530/910/920/930PG, ZM-920/930PG). The long battery life reduces running costs.

Measurement Parameters



ZS-910PG
ECG



ZS-920PG
ECG, Resp

* 8-lead ECG meas. available (ZM-920PG)



ZS-930PG
ECG, Resp, SpO₂

* LCD shows SpO₂
* 8-lead ECG meas. available (ZM-930PG)



ZS-940PG/ZM-940PG
NIBP, ECG, Resp, SpO₂

* LCD shows NIBP, SpO₂ and pulse rate

* NIBP measurement mode can be selected from manual, continuous and interval

* 8-lead ECG measurement available (ZM-940PG)



ZS-530PG
ECG, Resp, SpO₂

* Display shows ECG, heart rate, pulse rate, respiration rate, SpO₂

* ECG: one waveform from lead I, II or III

Related Products



Central monitor
CNS-9601



Telemetry system
WEP-5204/5208



Multiple patient receiver
ORG-9100/RU-960P

Options/Consumables

For full list of options/consumables, see the Technical Data separately available.

SpO₂ probe (BluPRO®)



Reusable type
P225H, finger, 0.6 m cord
P255F, finger, 1.6 m cord
P225G, multi-site



Disposable type
Adult to child, 4 types

This brochure may be revised or replaced by NIHON KOHDEN at any time without notice.

 **NIHON KOHDEN**

NIHON KOHDEN EUROPE GmbH
Raiffeisenstraße 10, 61191 Rosbach, Germany
Phone: +49 (0) 60 03 / 8 27-0, Fax: +49 (0) 60 03 / 8 27-5 99
Internet: www.nihonkohden.net, E-Mail: info@nke.de